



Extracting Forceps English Pattern




Extracting Forceps English Pattern




01

Extracting Forceps English Pattern




Extracting Forceps English Pattern



Extracting Forceps English Pattern



786-130
Fig. 51A
Upper roots
narrow beak



786-131
Fig. 51C
Upper roots



786-132
Fig. 52
Upper roots & canines
with wider beak



786-133
Fig. 52A
Upper teeth &
resection



786-134
Fig. 55
Separating upper
molars roots



786-135
Fig. 56
Separating lower
molars roots



786-136
Fig. 67
Upper wisdoms



786-137
Fig. 67A
Upper wisdoms

Extracting Forceps English Pattern



786-138
Fig. 68
Lower roots



786-139
Fig. 73
Lower molars
(Hawk's bill)



786-140
Fig. 74
Lower incisors
& roots



786-141
Fig. 74M
Lower incisors
& roots



786-142
Fig. 74N
Lower incisors
& roots



786-143
Fig. 75
Lower premolars
& canines



786-144
Fig. 76
Upper roots



786-145
Fig. 76N
Upper roots
(narrow beak)

02

Extracting Forceps English Pattern



786-146
Fig. 79
Lower wisdoms



786-147
Fig. 79A
Lower wisdoms



786-148
Fig. 86A
Lower molars



786-149
Fig. 86B
Lower molars



786-150
Fig. 86C
Lower molars


Extracting Forceps English Pattern



786-154
Fig. 94
Right upper molars



786-155
Fig. 95
Left upper molars



786-156
Fig. 99
Lower molars
(broken-down)



786-157
Fig. 99 1/4
Lower molars
(broken-down)



786-151
Fig. 87
Cowhorn beak




786-152
Fig. 89
Right upper molars



786-152/1
Fig. 22s
lower molars



786-153
Fig. 90
Left upper molars



786-158
Fig. 121
Upper wisdoms



786-159
Fig. 136
Upper centrals
& roots




786-160
Fig. 137
Lower centrals
& roots



786-161
Fig. 167
Universal pincer
for lower teeth
molars & premolars

Extracting Forceps English Pattern



786-162
Fig.168
Upper canines

786-163
Fig.169
Lower canines
& premolars

786-164
Fig.233
Lower roots
(parallel beaks)

786-165
Fig.235
Upper roots
(parallel beaks)



786-166
Fig. MD 1
Universal forceps for upper
wisdoms, premolars & roots

786-167
Fig. MD 2
Upper molars

786-168
Fig. MD 3
lower incisors, premolars
& roots

786-169
Fig. MD 4
Lowers molars

Extracting Forceps Children's Pattern



786-170

Fig. 158


Children's upper molars universal



786-171

Fig. 159

Children's upper teeth universal



786-172

Fig. 160


Children's lower molars universal(hawk's bill)



786-173

Fig. 161

Children's lower centrals & roots



786-174

Fig. 162

Children's upper centrals




786-175

Fig. 163

Children's upper centrals

03

Extracting Forceps English Pattern



786-176

Fig. 37


Children's upper incisors & canines



786-177

Fig. 39


Children's upper molars universal



786-178

Fig. 123


Children's lower teeth universal



786-179

Fig. 138

Children's upper universal & roots




786-180

Fig. 157

Children's upper molars universal

Wisdom Teeth Extracting Forceps Root Fragment Forceps



786-181

universal pattern for
upper teeth

786-182

universal pattern for
lower teeth

786-183

Fig. 22 1/2 L


Left lower wisdoms
& molars

786-184


Fig. 22 1/2 R

Right lower wisdoms
& molars

Tooth Extracting Forceps for children's with Spring




Tooth Extracting Forceps American Pattern




04


Tooth Extracting Forceps American Pattern



Tooth Extracting Forceps American Pattern




Tooth Extracting Forceps American Pattern



786-210
Fig. 88R
Right upper molars

786-211
Fig. 88L
Left upper molars

786-212
Fig. 99C
upper premolars,
molar & incisors




786-213
Fig. 101
Universal forceps
for premolars.

786-214
Fig. 150
Upper premolars,
incisors & roots

786-215
Fig. 150A
Upper premolars,
incisors & roots

786-216
Fig. 150S
Upper premolars,
incisors & roots
for children


Tooth Extracting Forceps American Pattern



786-217
Fig. 151
Lower premolars,
incisors & roots

786-218
Fig. 151A
Lower premolars,
incisors & roots

786-219
Fig. 151S
Lower premolars,
incisors & roots
for children



786-220
Fig. 210H
Upper wisdoms
with hook handle

786-221
Fig. 210S
Upper wisdoms
with straight handle

786-222
Fig. 222
Lower wisdoms

786-223
Fig. 286
Upper roots,
premolars
& incisors

05

Root Splinter Forceps Extraction Tweezer




786-224
Fig. 1
Lower roots

786-225
Fig. 2
Upper roots

786-226
Albrecht
16 cm


786-227
Albrecht
16 cm



786-228
Fig. 1
Root Fragments

786-229
Fig. 2
Root Fragments

Bone Rongeurs




786-230
Blumenthal
15 cm

786-231
Blumenthal 45°
14,5 cm

786-232
Blumenthal 90°
14 cm

786-233
Mead 1A
16,5 cm



786-234
Beyer
Multiple
Action
18 cm

786-235
Ruskin-mini
15,5 cm


786-236
Luer
15 cm

786-237
Friedman
14 cm

786-238
Mini-Friedman
12 cm


Root Elevators

Root Elevators



786-240 Fig.2B Bein
786-241 Fig.3B Bein
786-242 Fig.4B Bein
786-243 Fig.5B Bein
786-244 Fig.6B Bein
786-245 Fig.4SA Schmeckebier-Apexo
786-246 Fig.5SA Schmeckebier-Apexo

786-239 Fig.1B Bein
786-247 Fig.34S Seldin
786-248 Fig.34 Seldin
786-249 Fig.2S Seldin
786-250 Fig.3S Seldin
786-251 Fig.1F Flohr
786-252 Fig.2F Flohr
786-253 Fig.3F Flohr




786-255 Fig.72 Miller-Apexo
786-256 Fig.73 Miller-Apexo
786-257 Fig.74 Miller-Apexo
786-258 Fig.19
786-259 Fig.20
786-260 Fig.51 Curtis
786-261 Fig.52 Curtis

786-254 Fig.71 Miller-Apexo
786-262 Fig.1L Seldin
786-263 Fig.1R Seldin
786-264 Fig.4L Seldin
786-265 Fig.4R Seldin
786-266 Fig.11 Stout
786-267 Fig.11A Stout
786-268 Fig.44C Stout




Root Elevators

Root Elevators



786-270 Fig.B Cogswell
786-271 Fig.301 Apexo
786-272 Fig.302 Apexo
786-273 Fig.303 Apexo
786-274 Fig.301W Seldin
786-275 Fig.304W Seldin


786-269 Fig.A Cogswell
786-276 Fig.1H Heidbrink
786-277 Fig.2H Heidbrink
786-278 Fig.40H Heidbrink
786-279 Fig.41H Heidbrink
786-280 Fig.46
786-281 Fig.46R




786-283 Fig.45 Cryer
786-284 Fig.25 Cryer
786-285 Fig.26 Cryer
786-286 Fig.31 Cryer
786-287 Fig.32 Cryer
786-288 Fig.80
786-289 Fig.81

786-282 Fig.44 Cryer
786-290 Fig.81A
786-291 Fig.81E
786-292 Fig.190
786-293 Fig.191
786-294 Fig.23
786-295 Fig.24
786-296 Fig.98

Root Elevators




Root Elevators




07

Root Elevators




Pliers For Orthodontics & Prosthetics




Pliers For Orthodontics & Prosthetics




Pliers For Orthodontics & Prosthetics



Pliers For Orthodontics & Prosthetics



Orthodontics Instruments





786-001
De Bonding Plier
 Bracket / Adhesive remover with replaceable tool steel tip. Wedge between both edges of the base of the bracket and tooth surface. Easily lifting off the bracket or Bonding material. 49.50



786-002
Posterior Band Remover Long
 The padded tip rest on the Occlusal of the molar and the carbide tipped side removes the bracket or Bonding material. 49.50



786-003
Posterior Band Remover Small
 49.50



786-004
Tweed Loop Forming Plier
 3 Size carbide tip allows to make 3 size loops. Can hold wire firmly upto .022" x .028". The other side serrated tip gives none slip grip. 49.50

All Pliers Also Available In Box Joint

Orthodontic Instruments



786-005
How Plier
 Firmly hold the bracket tie wing for seating of the band. Carbide serrated tips for non slip grip. 49.50



786-006
Jara Back Plier
 A Plier for precise bending of wires. Precision grooves ensure accurate bending and closing of loops. Serrated tip for up to .020" wires. 49.50



786-007
Hard Wire Cutter
 Cuts hard wire up to .020" carbide insert for longer life. 49.50



786-008
Pin & Ligature Cutter
 Standard size lips designed to cut soft wire up to .051". 49.50

All Pliers Also Available In Box Joint



786-009
Angled Direct Bond Bracket Remover
 60° Angled tip allows for easy access to difficult anterior and posterior areas. 49.50



786-010
Weingart Utility Plier
 Matching deep serrated pads holds wire firmly at a convex working angle. Tapered narrow lip fits easily between brackets. 49.50

Orthodontic Instruments



786-011
Three Prong Clasps Adjusting Plier
 Rounded jaws close and adjust clasps without niking wire upto .030". 49.50



786-012
Light Wire Plier
 Wire forming Plier upto .061" the longer tapered beaks make it easier to bend small diameter loops. 3Grooves on square tip. 49.50

All Pliers Also Available In Box Joint



786-013
Bird Beak Plier
 1/2" Long tip carbide inserted round tip and flat ground tip gives non slip grip for wire up to .030". 49.50



786-014
Safety hold, small Distal end Cutter
 Narrow and light. 49.50

Orthodontic Instruments



786-015
Safety hold, Distal end Cutter
 Carbide tips for longer life and sharpness cuts and hold wire up to .022" x .028". 49.50





786-016
Mini Pin & Ligature Cutter
 Narrow tips cuts pins, Ligatures and elastic smoothly and soft wire up to .015". 49.50

All Pliers Also Available In Box Joint



Orthodontics Instruments

Gum Scissors




Gum Scissors


Gum Scissors




Cotton and Dressing Pliers




Needle Holders




Designed for Frontal & Palate Suture



Needle Holders



Towel Clamps, Sterilizing Forceps for Sterile Instruments



Elevators



786-460
Williger
16,5 cm



786-461
Williger
17 cm



786-462
König
15,5 cm



786-463
König
15,5 cm




786-464
Freer
19 cm



786-465
Freer
18,5 cm



786-466
Freer-Yasargin
18,5 cm



786-467
Kendzia
17 cm

Elevators




786-468
Söldin
19,5 cm




786-469
Söldin
19,5 cm



786-470
Chompret
15 cm



786-471 786-472



786-473 786-474



786-475
Freer
16 cm



786-476
Freer
16 cm



786-477
Partsch, mod
14,5 cm




786-478
Mead
16,5 cm



786-479
Mead
18 cm

12

Raspatories



786-480 3 mm
786-481 6 mm
17 cm



786-482 8 mm
Trelat
14 cm



786-483 5 mm
786-484 4.2 mm
König
15,5 cm
Williger
14 cm




786-485
Williger
16 cm




786-487 786-488
6 mm 8 mm

786-486 6 mm
16,5 cm



786-489
Obwegeser
17 cm



786-491 5 mm
786-492 8 mm

786-490 5 mm
20 cm



786-493
Tessier
22 cm

Chisels & Gouges



786-495
Partsch
17 cm



786-496
Freer
16 cm



786-497
Freer
16 cm



786-498
Lucas
16 cm



786-499
Lucas
16 cm



786-494
Lambotte
12,5 cm



786-500
Partsch
14 cm



786-501
Partsch
14 cm




786-502
Partsch
17 cm




786-503
Partsch
14 cm

MOUTH MIRRORS AND MIRRORS HANDLES


Cone Socket
786-507
 Regular, plain, 3, 4 and 5
786-508
 Regular, Magnifying, 3, 4 and 5
786-509
 Front-Surface, 3, 4 and 5



Mouth Mirrors Handles For Mouth Mirrors




Handles for Mouth Mirrors




13

Handles Reference Guide




Note:
 If you have preference different from the standard handle shown, please specify the optional handle number when ordering. If the tip/handle combination you desire is not shown, we will be glad to provide you with a quote for its preparation.

Explorers Endodontic Instruments




Periodontal pocket probes Endodontic Instruments




14

Explorers Probes





Explorers Endodontic Instruments



Filling Instruments



Filling Instruments




15

Filling Instruments


Filling Instruments Plastic Filling Instruments



Plastic Filling Instruments Amalgam Carvers




Amalgam Carvers Operative Instruments




16


Excavators




Cavity Preparation Instruments



Excavators




Excavators English Pattern




17

Bone Files



Bone Files




Periodontia Instruments




786-765 Fig. M9 Molt
 786-766 Fig.1 Hopkins
 786-767 Fig.2 Hopkins
 786-768 Fig.3 Mead
 786-769 Fig. PR3 Prichard
 786-770 Fig.15 Freer
 786-771 Fig.16 Freer

Bone Curettes Periodontia Instruments



786-772 Fig.1 Hemingway
 786-773 Fig.2
 786-774 Fig.3
 786-775 Fig.4
 786-776 Fig.85
 786-777 Fig.86
 786-778 Fig.87
 786-779 Fig.88
 Lucas
 Hemingway
 786-780 Fig.1S/2S Sugarman
 786-781 Fig.3S/4S Sugarman
 786-782 Fig.9/10 Schluger
 786-783 Fig.11/12 Buck
 786-784 Fig.3 Ochsenein
 786-785 Fig.4 Ochsenein
 786-786 Fig.36/37 Rhodes Back-Action


Periodontia Instruments



786-787 Fig.7 Goldman-Fox
 786-788 Fig.8 Goldman-Fox
 786-789 Fig.9 Goldman-Fox
 786-790 Fig.11 Goldman-Fox
 786-791 Fig.3/4 Buck
 786-792 Fig.5/6 Buck
 786-793 Fig.1/2 (15k/16k) Kirkland
 786-794 Fig.1/2 Orabn
 786-795 Fig.1/2 Sanders
 786-796 Fig.4/5 Sanders
 786-797 Fig.19/20 (University of South California) USC Towner
 786-798 Fig. CK3 Crane-Kaplan


Periodontia Instruments

curettes and scalers




786-800 Fig.2R/2L Columbia
 786-801 Fig.4R/4L Columbia
 786-802 Fig.13/14 Columbia
 786-803 Fig.411/412
 786-804 Fig.413/414
 786-805 Fig.1/2 Langer
 786-806 Fig.3/4 Langer
 786-799 Fig. HU 1/2 Curette
 786-807 Fig.1/2 Barnhart
 786-808 Fig.5/6 Barnhart
 786-809 Fig. ME-1 [10/11]
 786-810 Fig.3/4 Rule [23/24]
 786-811 Fig.11/14 Gracey Mesial-Distal
 786-812 Fig.12/13 Gracey Mesial-Distal
 786-813 Fig.13/14 Indiana University Anterior

Periodontia Instruments Gracey Curettes




Periodontia Instruments Scalers




19

Periodontia Instruments Scalers



Periodontia Instruments



CURETTES

ELEVATORS



786-876 786-877 786-878 786-879 786-880 786-881 786-882




786-883 Freer 786-884 West 786-885 786-886 (S23) Seldin 786-887 Ms.Molt 786-888 (6X) 786-889 (7X)

20

EXCAVATORS

EXCAVATORS





786-890 (20-L) 786-891 (30-L) 786-892 (3CHK-L) 786-893 (50-L) 786-894 (G160-L) 786-895 (3AS-L) 786-896 (6S-L) 786-897 (17S-L) 786-898 (238-L)



786-899 786-900 786-901 786-902 786-903 786-904 786-905 786-906 786-907 786-908 786-909



PROBE & EXPLORERS

SCALERS




Finishing Instruments


FILLING INSTRUMENTS




Wax & Modelling Instruments




Nerve Instruments




Alveotomy Files, Periosteal Elevators




Apical Fragment Ejectors



Bone Currettes





Periosteal Elevators, Bone Currettes




23


Scalpels - Gum Knives



Chisels and Gouges, Mallets




Wax Knives



786-982 Fig. 1 Fahnenstock cm. 12,5
 786-983 Fig. 1 Lessmann cm. 12,5
 786-984 Fig. 1 Con Scodellino cm. 12,5
 786-985 Fig. 2 Fahnenstock cm. 17
 786-986 Fig. 2 Lessmann cm. 17
 786-987 Fig. 2 Con Scodellino cm. 17

Spatulas for Plaster and Alginate



786-988
cm. 17

786-989
cm. 17 Gritman

786-990
cm. 19,5

786-991
cm. 21,5

786-992
cm. 20,5


786-993
cm. 19,5

786-994
cm. 20

786-995
cm. 22,5

24

Cement Spatulas Heidemann Separating Spatulas



786-996

786-997

786-998

786-999

786-1000

786-1001

786-1002

786-1003

786-1004
Fig.1

786-1005
Fig.2

786-1006
Fig.3

Wax Knives



786-1007
Fig. 6R

786-1008
Fig. 7R

786-1009


786-1010

786-1011

786-1012
Gritman


786-1013
Sachs

Wax and Modelling Carvers



786-1014 Fig.1 P.K. Thomas
 786-1015 Fig.2 P.K. Thomas
 786-1016 Fig.3 P.K. Thomas
 786-1017 Fig.4 P.K. Thomas
 786-1018 Fig.5 P.K. Thomas

Wax and Modelling Carvers



786-1019 zahle
 786-1020 zahle
 786-1021 zahle
 786-1022 Hylin
 786-1023 Hylin-Cermaica
 786-1024 Hylin
 786-1025 Hylin-Le Cron
 786-1026


25

Wax and Modelling Carvers



786-1027 Le Cron
 786-1028 Le Cron
 786-1029 Vehe
 786-1030 Beale
 786-1031 Roach
 786-1032 Le Cron
 786-1033 Le Cron
 786-1034 Gritman

Amalgam Carriers



02 mm 02.5 mm 03 mm
 786-1035 Small
 786-1036 Medium
 786-1037 Large
 786-1038 Jumbo
 786-1039 Amalgam well
 02 mm 02 mm 02 mm
 786-1040 Medium /Large
 786-1041 Medium /Small
 786-1042 Medium /Jumbo
 786-1043
 786-1044
 786-1045

Plaster Cutting Plier



786-1046

786-1047

786-1048

Instruments for Crown Work




786-1054
Pliers for measuring
crowns

786-1055
Steiner
Band Driver

786-1056
Crown and Band
Slitters

Crown Instruments Plaster Cutting Plier



786-1049
15 cm

Articulating
paper forceps


786-1050
Miller
15 cm

786-1051
Crown Spreader

786-1052
17 cm
786-1053
20 cm
Plaster Cutting Plier

26

Chip Blowers, Water Syringes Napkin Holder, Film Holder



786-1058

786-1059

786-1060

786-1057
Chip Blower

786-1064
Water Syringe

786-1065
Water Syringe

786-1061
Film Holder for
10 films
786-1062
6 films
786-1063
14 films

786-1066
Napkin holder with plastic band


786-1067
Napkin holder with metal chain

Matrix Retainers

786-1068 - 5mm
786-1069 - 6mm
786-1070 - 7mm
Ivory 8

786-1073
Siqueland-Narrow
786-1074
Siqueland-Wide


786-1076
Fig. 3
Tofflemire-Senior



786-1071
Ivory 8 N
786-1072
Ivory 8

786-1075
Fig. 1
Tofflemire-Universal

786-1077
Fig. 2
Tofflemire-Junior




786-1078
Fig. 2
Didner 8P

786-1079
Ivory

786-1080
Ivory

786-1081
Ivory
Separator

Crown Instruments



Straight, saw edge

Curved, saw edge

786-1082
Fig.1
10,5 cm

786-1083
Fig.2
10,5 cm

786-1084
Universal
12 cm

786-1085
For metal
786-1086
For wax

Universal


786-1087
Ivanson

786-1088
Fig.3
Crown spreaders

786-1089
Boley
14 cm

27

Crown Removers



786-1090
Crown removers
Automatic


with 5 Points

786-1091
Crown removers

786-1092
Crown removers

786-1093
Crown removers

Dental Cartridge Syringes



786-1094
Intraligamentary Syringe
1,8 ml

786-1095
Cartridge Syringe
2,2 ml
1,8 ml

786-1096
Aspirated Syringe
1,8 ml


786-1097
Aspirated Syringe
1,8 ml

786-1098
Cartridge Syringe
1,8 ml

Available Thread:- Inches & Metric



Soldering Tweezer

Bone Cutting Forceps



Rubber Dam Pliers

Articulators




Stainless Steel Impression Trays Made of 18/8 Stainless Steel



All the above impression trays are available in these 5 sizes

UPPER	All the above impression trays are available in these 5 sizes					LOWER
	U1	U2	U3	U4	U5	
A	62	60	53	54	52	L1
B	83	75	73	69	68	L2
						L3
						L4
						L5

All Dimension are in mm



Perforated Partial Set of 10
786-010wa


Stainless Steel Impression Trays Made of 18/8 Stainless Steel



All the above impression trays are available in these 3 sizes

UPPER	All the above impression trays are available in these 3 sizes			LOWER
	U1	U2	U3	
A	67	64	57	L1
B	81	71	70	L2
				L3

All Dimension are in mm



Perforated Set of 6
786-024wa

Stainless Steel Impression Trays Made of 18/8 Stainless Steel



All the above impression trays are available in these 4 sizes

UPPER	All the above impression trays are available in these 4 sizes				LOWER
	XL	L	M	S	
A	70	67	64	60	XL
B	79	76	70	73	L
					M
					S

All Dimension are in mm



Perforated Full Denture Set of 8
786-018wa

29


Stainless Steel Impression Trays Made of 18/8 Stainless Steel



All the above impression trays are available in these 4 sizes


UPPER	All the above impression trays are available in these 4 sizes				LOWER
	XL	L	M	S	
A	70	67	64	60	XL
B	79	76	70	73	L
					M
					S

All Dimension are in mm



Perforated Full Denture Set of 8
225-035wa

Stainless Steel Impression Trays Made of 18/8 Stainless Steel




All the above impression trays are available in these 4 sizes

Size	A	B	C
0	70	50	14
1	77	51	15
2	79	55	15
3	80	60	17


All Dimension are in mm

Stainless Steel Impression Trays Made of 18/8 Stainless Steel

Impression trays for crown & bridge work




Stainless Steel Impression Trays Made of 18/8 Stainless Steel



30

Rubber Dam Clamps



Mouth Gags



786-1212
Denhart
(12 cm)



786-1213 (16 cm)
786-1214 (20 cm)
Roser-Konig



786-1215
Heister
(13 cm)



786-1216
Molt
(11 cm)



786-1217
Molt
(14 cm)

Tongue Holding Forceps



786-1218
Young
(17 cm)



786-1219



786-1220
Collin
(17 cm)



786-1221
Guy
(19 cm)


Tongue Depressors, Cheek And Lip Retractors



786-1223
Hartmann
(15 cm)




786-1224
Tobold
(14.5 cm)



786-1225
Tobold
(14.5 cm)



786-1222
Andrew
(11 cm)




786-1226 to
1227
Wieder

Sternberg
786-1231
(14 cm)
786-1232
(16 cm)
786-1233
(18 cm)

786-1226
(14 cm)

786-1227
(15 cm)

Jennings
786-1234 (13 cm)
786-1235 (15 cm)
786-1236 (9 cm)
786-1237 (11 cm)



Sternberg
786-1228
(12 cm)
786-1229
(14 cm)
786-1230
(16 cm)



31

Retractors

786-1238 to 1239
Weitlaner-Wullstein
(13 cm)




1/1



786-1238



786-1239




1/1


786-1240
Mollison
(15 cm)




1/2




1/1



786-1241



786-1242




786-1244

786-1241 to 1244
Adson
(26.5 cm)


Holloware




Instruments Tray
With Lid
786-1245



Instruments Box
786-1246
8" x 12"
with lid




786-1247
Endodontic box
Also available in Solid pattern



Instruments Box
With Cover
786-1248
6" x 10"
without lid

Holloware



786-1249
Scalars Tray



Instruments Tray
with lid
786-1150
6" x 8"



Instruments Trays
with Or Without Lid

786-1151
6" x 8"
8" x 10"
10" x 12"

Holloware




Lotion Bowl
786-1252
2.25" x 2"
4.75" x 2"
5.5" x 2.5"
6.25" x 2.75"
7" x 3.5"

786-1253
80x35 mm
125x60 mm
Medicine Cup
Iodine Cup



786-1254
Cotton holder
60x50



786-1255
Forcep Jar
25 x 80 mm
30 x 100 mm
100 x 200 mm
75 x 173 mm

Holloware



Kidney tray
Traditional

786-1256
12"-10" - 8" - 6"



Deep

786-1257
6" x 9" - 8" x 12"



Shallow

786-1258
10" x 10"



786-1259
Spirit lamp
50cc,100cc



786-1260
Spirit lamp
50cc,100cc